


# X-ray direct-detection CCDs

#### Preliminary Short Form Datasheet



CCD3200A Direct Detection CCD Si-CCD chip on sensor-board

### **General Description:**

Direct Detection CCDs convert x-rays into electrical carriers directly in the Silicon, without using a scintillator to transform the x-ray radiation into visible light first.

#### X-ray energy bands:

0.12 eV – 12 keV (soft x-rays): Several local detection maxima with high x-ray detection efficiency up to > 80 %, e.g. at 12 keV (without window)

12 keV – 120 keV (hard x-rays): can also be well detected with a scintillator in front of the CCD.

#### **Design Examples:**

### Front Illuminated (FI):

- Highest quantum efficiency for E > 4 KeV
- Number of Pixels: 4k x 4k
- Pixel Size: 24μm x 24μm
- No. of outputs: 16
- Frame Rate: 2 fps
- X-ray resolution: 22µm (PSF at FWHM)
- Lower Cost than BI

### Back Illuminated (BI):


- Highest quantum efficiency for E < 12 KeV
- Number of Pixels: 4k x 4k
- Pixel Size: 24μm x 24μm
- No. of outputs: 16
- Frame Rate: 2 fps
- X-ray resolution: e.g. 20 (40)µm at 20 (10) keV (PSF at FWHM)

# X-ray specific properties:

- Excellent X-ray detection efficiency combined with highest spatial resolution
- Detection efficiency up to > 80 %, e.g. at 12 keV (without window)
- Spatial Resolution up to 4k x 4k (16 MPixels) monolithic,
- X-ray resolution down to 50μm and better (20 μm) (FWHM of PSF)
- Low Radiation Damage by adequate CCDdesign

# **Si-CCD-specific properties:**

- 100% Fill Factor!
- Deep depletion high resistivity silicon
  > 10 k Ωcm, e.g. 700 µm thick
- Back illuminated or Front illuminated
- Pixel Sizes: 7.5 μm 150 μm (24 μm x 24 μm at 4k x 4k resolution)
- Chip-Dimensions: up to 96 mm x 96 mm (monolithic)
- Pixel-Resolution up to 10 k x 10k (111 Mpixel) (monolithic)
- Low Crosstalk between channels



CCD3200A-Wafer: Chip-Dimension 96mm x 96mm

# **Applications:**

- Synchrotron Research and Measurements
- Crystallography
- X-ray astronomy
- X-ray microscopy and spectroscopy
- High Energy Physics
- Medical research

Copyright © 2014 ANDANTA GmbH; Version 0.2. dtd. February 8, 2014. All rights reserved. The information contained in this document has been summarized to the best of our knowledge. However, no responsibility is accepted for the consequences of any use thereof. Furthermore, the information provided may be changed without explicit notice.

ANDANTA GmbH Detektortechnologie Ilzweg 7+9 • 82140 Olching/Deutschland Tel:+49814241058-0 • Fax:+49814241058-29 e-mail: epost@andanta.de • www.andanta.de